FPSA

Free Prostate Specific Antigen (CLIA)

Order Information

Catalog No.	Package Size
FPSA111	2×50 tests
FPSA112	2×100 tests

Intended Use

The CL-series FPSA assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative determination of free prostate specific antigen (Free PSA, FPSA) in human serum.

The CL-series FPSA assay can be used with the t-PSA as an aid in distinguishing prostate cancer from benign prostatic conditions. Prostate biopsy is required for the diagnosis of prostate cancer.

Summary

Prostate cancer is one of the most common cancers in males. Early detection is very important for medical treatment of prostate cancer patients, and it requires a simple, safe, and inexpensive test for the disease.^{1, 2, 3}

Prostate Specific Antigen (PSA) is a single chain glycoprotein produced in the glandular epithelium of the prostate. It is secreted into the seminal fluid in high concentration.² PSA is a serine protease with chymotrypsin-like activity. A major function of PSA is the proteolytic cleavage of gel-forming proteins I and II in the seminal fluid, resulting in the liquefaction of the seminal gel and increased sperm mobility. Normally PSA concentration in spermiatic fluid is high, but its concentration in serum is very low.⁴ Increasing levels of serum PSA are associated with pathologic conditions of prostate (e.g. Prostatitis, Benign Prostatic Hyperplasia or Prostatic Cancer).

PSA exists in three major forms in blood: free PSA, PSA complexed with the serine protease inhibitor, alpha-1-antichymotrypsin (PSA-ACT), and PSA complexed with alpha-2-macroglobulin that is lack of immune reaction.^{5, 6, 7} Only free PSA and PSA-ACT are detectable in current immunoassay, and they are defined together as total PSA.⁸ Generally, the lower the ratio of serum free PSA/total PSA, the higher the possibility of prostate diseases in men.^{7, 8}

Assay Principle

The CL-series FPSA assay is a two-site sandwich assay to determine the level of free PSA.

In the first step, sample, paramagnetic microparticle coated with monoclonal anti-free PSA antibody (mouse) and monoclonal anti-PSA antibody (mouse)-alkaline phosphatase conjugate are added into a reaction vessel. After incubation, free PSA present in the sample binds to both anti-PSA antibody coated microparticle and anti-PSA antibody alkaline phosphatase-labeled conjugate to form a sandwich complex. Microparticle is magnetically captured while other


unbound substances are removed by washing.

In the second step, the substrate solution is added to the reaction vessel. It is catalyzed by anti-PSA antibody (mouse)-alkaline phosphatase conjugate in the immunocomplex retained on the microparticle. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier built into the system. The amount of free PSA present in the sample is proportional to the relative light units (RLUs) generated during the reaction. The Free PSA concentration can be determined via a calibration curve.

Reagent Components

Ra	Paramagnetic microparticles coated with monoclonal anti-free PSA antibody (mouse) in TRIS buffer with preservative.
Rb	Monoclonal anti-PSA antibody (mouse)-alkaline phosphatase conjugate in PBS buffer with preservative.

The position of each reagent component is shown in the figure below (front view on the left and top view on the right):

Storage and Stability

The unopened FPSA (CLIA) reagent kit is stable up to the stated expiration date when stored at 2-8°C. The FPSA (CLIA) reagent kit can be stored onboard and used for a maximum of 28 days after opening at 2-8°C.

Reagent Preparation

Ra: Ready to use; Rb: Ready to use

Materials Required but not Provided

Mindray CL-series Chemiluminescence Immunoassay Analyzer.

Cat.No.FPSA211: Mindray Free PSA Calibrators, 1×2.0 mL for each of calibrator C0, C1 and C2.

Cat.No.TML311/TML312/TML313/TML314: Mindray Tumor Marker Multi Control (L), 6×5.0 mL/12×5.0 mL/1×5.0 mL/3×5.0 mL.

Cat.No.TMH311/TMH312/TMH313/TMH314: Mindray Tumor Marker Multi Control (H), 6×5.0 mL/12×5.0 mL/1×5.0 mL/3×5.0 mL.

Cat.No.TML321/TML322/TML323/TML324: Mindray Tumor Marker Multi Control (L), 1×5.0 mL/3×5.0 mL/6×5.0 mL/12×5.0 mL.

Cat.No.TMH321/TMH322/TMH323/TMH324: Mindray Tumor Marker Multi Control (H), 1×5.0 mL/3×5.0 mL/6×5.0 mL/12×5.0 mL.

Cat.No.WB411: Mindray Wash Buffer, 1×10 L.

Cat.No.CS511/CS512: Mindray Substrate Solution, 4×115 mL/4×75 mL.

Mindray Reaction Vessels.

Applicable Instrument

Mindray CL-series Chemiluminescence Immunoassay Analyzer.

Specimen Collection and Preparation

Human serum is recommended for this assay.

Centrifuge the specimens after clot formation is complete. Transfer the supernatants into tubes for storage or test within two hours after centrifugation.

Specimens should be tested as soon as possible after sample collection. If testing is not completed within 8 hours, specimens should be tightly capped and refrigerated at 2-8°C. If testing will be delayed for more than 72 hours, specimens should be frozen at -20°C or below.

Avoid repeated freeze and thaw cycles.

Assay Procedure

For optimal performance of this assay, operators should read the related system operation manual carefully, to get sufficient information such as operation instructions, sample preservation and management, safety precaution, and maintenance. Prepare all required materials for the assay as well.

Before loading the FPSA (CLIA) reagent kit on the machine for the first time, unopened reagent bottle should be inverted gently for at least 30 times to resuspend the microparticles that have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been resuspended. If the microparticles remain adhered to the bottle, continue inverting until the microparticles have been completely resuspended. If the microparticles cannot be resuspended, it is recommended not to use this bottle of reagent. Contact Mindray Customer Service for help. Do not invert opened reagent bottle.

This assay requires 10 µL of sample for a single test. This volume does not include the dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the system operation manual and specific requirement of the assay to determine the minimum sample volume.

Calibration

CL-series FPSA (CLIA) has been standardized against the WHO International Standard Prostate-Specific Antigen Free (NIBSC code: 96/668).

The specific information of master calibration curve of FPSA (CLIA) reagent kit is stored in the two-dimensional barcode attached in the reagent pack. It's used together with calibrators for the calibration of the specific reagent lot. When performing the calibration, scan the information of master calibration curve from the barcode into the system first, and then use the calibrators at three levels. Valid calibration curve is required before any FPSA test. Recalibration is recommended every 4 weeks, or when a new reagent lot is used, or the

quality controls are out of specified range. For detailed instruction of calibration, refer to the system operation manual.

Quality Control

It is recommended that quality controls should be run once every 24 hours if the tests are in use, or after every calibration. The quality control frequency should be adapted to each laboratory's individual requirements. The recommended two levels of quality controls for this assay are Mindray Tumor Marker Multi Control (L) and Tumor Marker Multi Control (H).

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Examine the assay system referring to the system operation manual. If the quality control results are still out of the specified range, please contact Mindray Customer Service for help.

Calculation

The analyzer automatically calculates the analyte concentration of each sample on the master calibration curve read from the barcode, and a 4-Parameter Logistic Curve Fitting (4PLC) with the relative light units (RLUs) generated from three level calibrators of defined concentration values. The results are shown in the unit of ng/mL.

Conversion factors: ng/mL × 1 = µg/L

Expected values

An extensive study on a cohort of 304 healthy males has determined the reference range of CL-series FPSA assay.

Category	N	97.5% Percentile
Male	304	1.0 ng/mL

Due to the variation in geography, race, sex, and age, it is highly recommended that each laboratory should establish its own reference range.

Limitation

The upper limit of this assay is 30 ng/mL. A specimen with free PSA concentration lower than the upper limit can be quantitatively determined, while specimen with a concentration higher than the upper limit will be reported as >30 ng/mL.

The concentration of free PSA in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity. The assay results should be used in conjunction with other data, such as symptoms, results of other tests, clinical history, etc.

Specimen from individuals who have been exposed to mouse monoclonal antibodies may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values with assay kits employing mouse monoclonal antibodies. However, no obvious interference of HAMA was observed in this assay.

Performance Characteristics**Analytical Sensitivity/Limit of Detection**

The FP5A (CLIA) reagent kit has an analytical sensitivity of ≤ 0.01 ng/mL. Analytical sensitivity is defined as the lowest concentration of analyte that can be differentiated from a sample that contains no analyte. It is defined as the Free PSA concentration at two standard deviations above the mean RLU from 20 measurements of an analyte-free sample.

Reportable Range

Reportable range is defined by the analytical sensitivity and the upper limit of the master calibration curve. The reportable range of FP5A (CLIA) reagent kit is 0.01-30 ng/mL.

Specificity

Hemoglobin up to 500 mg/dL, bilirubin up to 20 mg/dL, triglycerides up to 1500 mg/dL, and total protein up to 10.0 g/dL will not interfere with the CL-series FP5A assay. These substances show less than 10% interferences at indicated concentration.

No obvious interference was observed from rheumatoid factor up to 400 IU/mL or antinuclear antibody up to 2000 U/L.

Interference substance tests were performed on 10 commonly used pharmaceuticals. These compounds showed less than 10% interference in the CL-series FP5A assay at the levels indicated below.

Tested Compound	Interference Substance Concentration
Leuporelin acetate	100 µg/mL
Cyclophosphamide	700 µg/mL
Finasteride	370 ng/mL
Megestrol	2.4 mg/dL
Methoposterin	30 µg/mL
Flutamide	10 µg/mL
Doxorubicin hydrochloride	16 µg/mL
Aspirin	0.5 mg/mL
Biotin	50 ng/mL
Diethyl oestrodione	2 µg/mL

Mindray Free PSA Calibrator C0 were supplemented with other tumor markers, such as alpha-fetoprotein (AFP), cancer antigen 125 (CA125), cancer antigen 15-3 (CA15-3), carbohydrate antigen 19-9 (CA19-9), carcinoembryonic antigen (CEA) and ferritin (FERR) at specific levels indicated in the table below. No obvious cross reactivity was observed as all the results were ≤ 0.5 ng/mL. The results are stated in the table below.

Tumor Marker	Cross-reactant Concentration	Reported FP5A (ng/mL)	Acceptance Criteria
AFP	1000 ng/mL	0.00	Reported FP5A
CA125	1000 U/mL	0.00	

CA15-3	100 U/mL	0.00	≤ 0.5 ng/mL
CA19-9	1000 U/mL	0.00	
CEA	1000 ng/mL	0.00	
FERR	1000 ng/mL	0.00	

High Dose Hook

For the CL-series FP5A assay, no high dose hook effect was observed when samples containing up to approximately 30,000 ng/mL of free PSA were assayed.

Representative performance data on CL-2000i analyzers are given below. For precision, performance data on CL-2000i and CL-1000i are both given.

Accuracy

A sample of WHO International Standard Prostate-Specific Antigen Free (NIBSC code: 96/668) with traceable and defined value was used to verify the accuracy of this assay. The results showed that the relative deviation was less than $\pm 10\%$. The results are listed in the following table.

Sample	Measured FP5A Value (ng/mL)	Defined FP5A Value (ng/mL)	Relative Deviation
WHO PSA FREE	10.37	10.12	2.50%

Precision

The CL-series FP5A assay is designed to have a precision of $\leq 10\%$ (within-device CV). Precision was determined by following National Committee for Clinical Laboratory Standards (NCCLS) Protocol EP5-A2. Two levels of quality controls were tested in duplicate in two separate runs per day for a total of 20 days, using a single lot of reagents and a single calibration curve. The precision data are summarized in the table below.

Performance data on CL-2000i				
Sample	Mean FP5A (ng/mL)	Within-run CV	Between-run CV	Within-Device CV
1	0.93	3.99%	2.89%	8.61%
2	11.46	1.28%	2.45%	5.53%

Performance data on CL-1000i				
Sample	Mean FP5A (ng/mL)	Within-run CV	Between-run CV	Within-Device CV
1	1.03	3.14%	1.88%	3.14%
2	10.54	3.26%	2.27%	3.26%

Linearity

A high concentration FP5A sample (approximately 20 ng/mL) was mixed with a low concentration sample (<0.01 ng/mL) at different ratios, generating a series of dilutions. The free PSA of each dilution was determined using the Mindray CL-Series FP5A Assay. Linearity was demonstrated in the range of 0.01

ng/mL to 20 ng/mL. The correlation coefficient r is ≥ 0.9900 . The linearity data are summarized in the table below.

Concentration (ng/mL)	1	2	3	4	5	6
Expected FP5A	0.0020	6.03	12.05	18.07	24.10	30.12
Measured FP5A	0.0020	6.28	12.79	18.69	24.18	30.12

Method Comparison

The Mindray CL-Series FP5A Assay was compared to a commercially available diagnostic kit in a correlation study with about 40 specimens. The statistical data obtained by Deming computing mode are shown in the table below.

Concentration Range (ng/mL)	Slope	Intercept	Correlation Coefficient
0.01-30	0.9218	0.1259	0.9973

Warnings and Precautions

- For in vitro diagnostic use only.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.


Warning: May cause an allergic skin

reaction.

Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/ protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of water. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.

- Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different manufacturers on Mindray system, or using Mindray reagent kits on other systems.
- Do not use reagent kits beyond the expiration date.
- Do not use reagents mixed from different reagent lots.
- Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
- Reagent pack opened for more than 28 days is not recommended for use.
- Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
- All the specimen and reaction wastes should be considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.
- The Material Safety Data Sheet (MSDS) is available upon request.

Graphical Symbols**References**

- Henttu P, Vihko P. Prostate-specific Antigen and Human Glandular Kallikrein: Two Kallikreins of the Human Prostate. *Ann Med*, 1994, 26(3):157-164.
- Prestigiacomo AF, Stamey TA. Clinical usefulness of free and complexed PSA. *Clin Lab Invest Suppl*, 1995, 221: 32-34.
- Christensson A, Laurell CB, Lilja H. Enzymatic activity of prostate-specific antigen and its reactions with extracellular serine protease inhibitors. *Eur J Biochem*, 1990, 194(3): 755-763.
- Armbruster DA. Prostate Specific Antigen: Biochemistry, Analytical Methods and Clinical Application. *Clin Chem*, 1993, 39(2):181-195.
- McCormack RT, Rittenhouse HG, Finlay JA, Sokoloff RL, Wang TJ, Wolfert RL, Lilja H, Oesterling JE. Molecular forms of prostate-specific antigen and the human kallikrein gene family: A new era. *Urology*, 1995, 45(5): 729-744.
- Zhou AM, Tewari PC, Bluestein BI, Caldwell GW, Larsen FL. Multiple forms of prostate specific antigen in serum: differences in immunorecognition by monoclonal and polyclonal assays. *Clin Chem*, 1993, 39(12):2483-2491.
- Zhang WM, Leinonen J, Kalkkinen N, Dowell B, Stenman UH. Purification and Characterization of Different Molecular Forms of Prostate-Specific Antigen in Human Seminal Fluid. *Clin Chem*, 1995, 41(11):1567-1573.
- Lilja H, Christensson A, Dahlen U, Matikainen M-T, Nilsson O, Petterson K, et al. Prostate Specific Antigen in Human Serum occurs predominantly in Complex with Alpha-1-Antichymotrypsin. *Clin Chem*, 1991, 37(9):1618-1625.

© 2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
All rights Reserved.

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Address: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

E-mail Address: service@mindray.com
Website: www.mindray.com

Tel: +86-755-26582888

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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- Henttu P, Vihko P. Prostate-specific Antigen and Human Glandular Kallikrein: Two Kallikreins of the Human Prostate. *Ann Med*, 1994, 26(3):157–164.
- Prestigiacomo AF, Stamey TA. Clinical usefulness of free and complexed PSA. *Clin Lab Invest Suppl*, 1995, 221: 32–34.
- Christensson A, Laurell CB, Lilja H. Enzymatic activity of prostate-specific antigen and its reactions with extracellular serine protease inhibitors. *Eur J Biochem*, 1990, 194(3): 755–763.
- Armbruster DA. Prostate Specific Antigen: Biochemistry, Analytical Methods and Clinical Application. *Clin Chem*, 1993, 39(2):181–195.
- McCormack RT, Rittenhouse HG, Finlay JA, Sokoloff RL, Wang TJ, Wolfert RL, Lilja H, Oesterling JE. Molecular forms of prostate-specific antigen and the human kallikrein gene family: A new era. *Urology*, 1995, 45(5): 729–744.
- Zhou AM, Tewari PC, Bluestein BI, Caldwell GW, Larsen FL. Multiple forms of prostate specific antigen in serum: differences in immunorecognition by monoclonal and polyclonal assays. *Clin Chem*, 1993, 39(12):2483–2491.
- Zhang WM, Leinonen J, Kalkkinen N, Dowell B, Stenman UH. Purification and Characterization of Different Molecular Forms of Prostate-Specific Antigen in Human Seminal Fluid. *Clin Chem*, 1995, 41(11):1567–1573.
- Lilja H, Christensson A, Dahlen U, Matikainen M-T, Nilsson O, Pettersson K, et al. Prostate Specific Antigen in Human Serum occurs predominantly in Complex with Alpha-1-Antichymotrypsin. *Clin Chem*, 1991, 37(9):1618–1625.

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
Símbolos gráficos

Referencias

1. Henttu P, Vihko P. Prostate-specific Antigen and Human Glandular Kallikrein: Two Kallikreins of the Human Prostate. Ann Med, 1994, 26(3):157-164.
2. Prestigiacomo AF, Stamey TA. Clinical usefulness of free and complexed PSA. Clin Lab Invest Suppl, 1995, 221: 32-34.
3. Christensson A, Laurell CB, Lilja H. Enzymatic activity of prostate-specific antigen and its reactions with extracellular serine protease inhibitors. Eur J Biochem, 1990, 194(3): 755-763.
4. Armbruster DA. Prostate Specific Antigen: Biochemistry, Analytical Methods and Clinical Application. Clin Chem, 1993, 39(2):181-195.
5. McCormack RT, Rittenhouse HG, Finlay JA, Sokoloff RL, Wang TJ, Wolfert RL, Lilja H, Oesterling JE. Molecular forms of prostate-specific antigen and the human kallikrein gene family: A new era. Urology, 1995, 45(5): 729-744.
6. Zhou AM, Tewari PC, Bluestein BI, Caldwell GW, Larsen FL. Multiple forms of prostate specific antigen in serum: differences in immunorecognition by monoclonal and polyclonal assays. Clin Chem, 1993, 39(12):2483-2491.
7. Zhang WM, Leinonen J, Kalkkinen N, Dowell B, Stenman UH. Purification and Characterization of Different Molecular Forms of Prostate-Specific Antigen in Human Seminal Fluid. Clin Chem, 1995, 41(11):1567-1573.
8. Lilja H, Christensson A, Dahlen U, Matikainen M-T, Nilsson O, Petterson K, et al. Prostate Specific Antigen in Human Serum occurs predominantly in Complex with Alpha-1-Antichymotrypsin. Clin Chem, 1991, 37(9):1618-1625.

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
Grafiksel Semboller

Referanslar

- Henttu P, Vihko P. Prostate-specific Antigen and Human Glandular Kallikrein: Two Kallikreins of the Human Prostate. Ann Med, 1994, 26(3):157-164.
- Prestigiacomo AF, Stamey TA. Clinical usefulness of free and complexed PSA. Clin Lab Invest Suppl, 1995, 221: 32-34.
- Christensson A, Laurell CB, Lilja H. Enzymatic activity of prostate-specific antigen and its reactions with extracellular serine protease inhibitors. Eur J Biochem, 1990, 194(3): 755-763.
- Armbuster DA. Prostate Specific Antigen: Biochemistry, Analytical Methods and Clinical Application. Clin Chem, 1993, 39(2):181-195.
- McCormack RT, Rittenhouse HG, Finlay JA, Sokoloff RL, Wang TJ, Wolfert RL, Lilja H, Oesterling JE. Molecular forms of prostate-specific antigen and the human kallikrein gene family: A new era. Urology, 1995, 45(5): 729-744.
- Zhou AM, Tewari PC, Bluestein BI, Caldwell GW, Larsen FL. Multiple forms of prostate specific antigen in serum: differences in immunorecognition by monoclonal and polyclonal assays. Clin Chem, 1993, 39(12):2483-2491.
- Zhang WM, Leinonen J, Kalkkinen N, Dowell B, Stenman UH. Purification and Characterization of Different Molecular Forms of Prostate-Specific Antigen in Human Seminal Fluid. Clin Chem, 1995, 41(11):1567-1573.
- Lilja H, Christensson A, Dahlen U, Matikainen M-T, Nilsson O, Petterson K, et al. Prostate Specific Antigen in Human Serum occurs predominantly in Complex with Alpha-1-Antichymotrypsin. Clin Chem, 1991, 37(9):1618-1625.

© 2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 Çin

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Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant EC : Shanghai International Holding Corp. GmbH (Europe)