CK

Generic Name : Creatine Kinase Kit（IFCC Method）
Abbreviated name : CK
Order Information

Cat. No.	Package size
CK0102	R1 2×40 mL + R2 1×20 mL
CK1102	R1 2×35 mL + R2 1×18 mL
CK2102	R1 1×15 mL + R2 1×8 mL
CK0103	R1 4×38 mL + R2 2×20 mL
CK0104	R1 3×45 mL + R2 3×12 mL
CK1104	R1 3×42 mL + R2 3×12 mL
CK0105	R1 2×250 mL + R2 1×125 mL

Intended use

In vitro test for the quantitative determination of creatine kinase (CK) activity in serum and plasma on photometric systems.

Summary [1]

Creatine kinase (CK) is an enzyme, which consists of isoenzymes mainly of the muscle (CK-M) and the brain (CK-B). CK exists in serum in dimeric form as CK-MM, CK-MB, and CK-BB and as macroenzyme. Elevated CK values are observed in cardiac muscle damages and in skeletal muscle diseases. Measurement of CK is used especially in conjunction with CK-MB for diagnosis and monitoring of myocardial infarction.

Method

International Federation of Clinical Chemistry and Laboratory Medicine （IFCC） method

Reaction Principle

$$\text{Creatine phosphate} + \text{ADP} \xrightleftharpoons{\text{CK}} \text{Creatine} + \text{ATP}$$

$$\text{ATP} + \text{Glucose} \xrightleftharpoons{\text{HK}} \text{ADP} + \text{Glucose-6-phosphate}$$

$$\text{Glucose-6-phosphate} + \text{NADP}^+ \xrightleftharpoons{\text{G-6-PDH}} \text{Gluconate-6-phosphate} + \text{NADPH} + \text{H}^+$$

Creatine kinase (CK) catalyzes the phosphorylation of ADP, in the presence of creatine phosphate, to form ATP and creatine. The catalytic concentration is determined from the rate of NADPH formation, measured at 340 nm, by means of the hexokinase (HK) and glucose-6-phosphate dehydrogenase (G-6-PDH) coupled reactions.

Reagents

Components and Concentrations

R1:	Imidazole buffer	100 mmol/L
	Glucose	20 mmol/L
	N-acetylcysteine(NAC)	0.2 mmol/L

	Magnesium acetate	10 mmol/L
	EDTA	2 mmol/L
	NADP	2 mmol/L
	AMP	5 mmol/L
	HK	>4 U/mL
	Creatine phosphate	30 mmol/L
R2:	ADP	2 mmol/L
	G-6-PDH	>2.8 U/mL

Warnings and Precautions

1.For in vitro diagnostic use.
2.Take the necessary precautions for the use of laboratory reagents.
3.Preservative contained. Do not swallow. Avoid contact with skin and mucous membranes.
4.Disposal of all waste material should be in accordance with local guidelines.
5.Material safety data sheet is available on request for professional users.

Reagent Preparation

R1 and R2 are ready to use.

Storage and stability

Stable up to expiry date indicated on the label, when stored unopened at 2-8℃ and protected from light.
Once opened, the reagents are stable for 30 days when refrigerated on the analyzer or refrigerator.
Contamination of the reagents must be avoided.
Do not freeze the reagents.

Reagent Blank Absorbance

The absorbance of reagent blank at 340 nm should be <0.5 A.

Materials required but not provided

1.Calibrator and controls as indicated below.
2.NaCl solution 9 g/L.
3.General laboratory equipments.

Specimen Collection and preparation [2]

1.Serum, EDTA plasma is suitable for samples. Whole blood, hemolysis and urine are not recommended for use as a sample. Freshly drawn serum is preferred specimen.
2.Use the suitable tubes or collection containers and follow the instruction of the manufacturer; avoid effect of the materials of the tubes or other collection containers.
3.Centrifuge samples containing precipitate before performing the assay.
4.Stability: 2 days at 20-25℃
1 week at 2-8℃

4 weeks at -20℃

Assay procedure

	Blank	Sample
Reagent 1	1000 µL	1000 µL
Dist. water	25 µL	—
Sample	—	25 µL
Mix, incubate for 3 min. at 37℃, then add:		
Reagent 2	250 µL	250 µL
Mix thoroughly, incubate at 37℃ for 3 min., and then read the absorbance change value within 1-3 min.		
ΔA/min = [ΔA/min sample]- [ΔA/min blank]		

Application sheets for BS series analyzers are available in this document. Please refer to the appropriate operation manual for the analyzer-specific assay instructions.

Calibration

1.It is recommended to use the Human multi-calibrator from Mindray and 9 g/L NaCl for two-point calibration. Traceability of the multi-calibrator can refer to the calibrator instructions for use of Mindray Company.

2.Calibration frequency:

After reagent lot changed.

As required following quality control procedures.

Quality control

At least two levels of control material should be analyzed with each batch of samples. In addition, these controls should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or troubleshooting procedures as detailed in the appropriate system manual.

We recommend using the Human Assayed Control made by Mindray to verify the performance of the measurement procedure; other suitable control material can be used in addition.

Each laboratory should establish its own internal quality control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

Calculation

The analyzer calculates the CK activity of each sample automatically with a specified valid calibration factor from calibration process.

Conversion factor of conventional units (U/L) into S.I. units (µkat/L):

1 U/L = 16.67×10^{-3} µkat/L

1 µkat/L = 60 U/L

Reference Intervals [3,4,5]

Each laboratory should establish its own reference intervals based upon its patient population. The reference intervals measured at 37℃ listed below

were taken from literature:

Sample Type		Conventional Units	S.I. Units
Serum / Plasma	Male	≤ 171 U/L	≤ 2.85 μkat/L
	Female	≤ 145 U/L	≤ 2.41 μkat/L

Performance Characteristics

Representative performance data obtained from Mindray system is given below. Results may vary if a different instrument, an individual laboratory or a manual procedure is used.

Interferences/Specificity

The following substances were tested for interference with this methodology. Criterion: Recovery within ±10 % of initial value.

Substance	Level Tested	Observed Effect
Ascorbic acid	30 mg/dL	NSI*
Bilirubin	40 mg/dL	NSI
Lipemia	500 mg/dL	NSI
Hemoglobin	200 mg/dL	NSI

* NSI: No Significant Interference (within ± 10 %)

Linearity range

The Mindray System (Mindray BS series analyzers / Mindray CK Reagent) provides the following linearity range

Sample Type	Conventional Units	S.I. Units
Serum / Plasma	5-1000 U/L	0.08-16.67 μkat/L

If the value of sample exceeds 1000 U/L, the sample should be diluted with 9 g/L NaCl solution (e.g. 1+ 9) and rerun; the result should be multiplied by 10.

Sensitivity/Detection Limit

The lowest measurable CK activity that can be distinguished from zero is 5 U/L (0.08 μkat/L) with 99.7% confidence.

Precision

Precision performance using the CLSI Approved Guideline EP5-A2 to assay serum control appears in the table below[6]. U:U/L

Type of Precision	Level Ⅱ			Level Ⅲ		
	Mean	SD	CV %	Mean	SD	CV %
Within-run	211.4	1.52	0.72	527.3	2.31	0.44
Between-run		1.71	0.81		3.27	0.62
Between-day		1.53	0.72		2.46	0.47
Within-device		2.75	1.30		4.70	0.89

Method Comparison

A comparison between Mindray system (Mindray BS series analyzers /Mindray CK reagent) (y) and Hitachi/Roche system (Hitachi /Roche CK) (x) using 40 samples gave following correlation (U/L): $y=1.0199x+5.9196$, R^2 =

0.9998 Details of the comparison experiments are available on request.

References

1.Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1999.
2.Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes;Preanalytical Variables. Brochure in: Samples: From the Patientto the Laboratory. Darmstadt: GIT Verlag 1996.
3.Schumann, G., et al., Clin. Chem. Lab. Med., 2002;40(6) pp 635 642.
4.National Committee for Clinical Laboratory Standards, How to Define, Determine, and Utilize Reference Intervals in the Clinical Laboratory, Approved Guideline, CLSI publication C28-A, Villanova, PA(1995).
5.Tietz, N. W., ed., Fundamentals of Clinical Chemistry, 3rd Edition, W.B. Saunders, Philadelphia, PA (1987).
6.CLSI. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP5-A2 [ISBN 1-56238-542-9. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087 USA, 2008.

Graphical symbols


In Vitro Diagnostic medical device

Batch Code

European Conformity

Authorized representative in the European Community

Use By

Consult Instructions for use

Temperature Limit

Manufacturer

Catalogue number

© 2015 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Address:Mindray building, Keji 12th Road South, Hi-tech Industrial Park,Nanshan, Shenzhen, 518057 P.R.China
E-mail Address: service@mindray.com
Website: www.mindray.com
Tel:+86-755-81888998 **Fax**:+86-755-26582680
EC-Representative: Shanghai International Holding Corp. GmbH(Europe)
Address: Eiffestraβe 80, Hamburg 20537, Germany
Tel: 0049-40-2513175 **Fax**:0049-40-255726

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Referências

1.Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1999.
2.Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the Patientto the Laboratory. Darmstadt: GIT Verlag 1996.
3.Schumann, G., et al., Clin. Chem. Lab. Med., 2002;40(6) pp 635 642.
4.National Committee for Clinical Laboratory Standards, How to Define, Determine, and Utilize Reference Intervals in the Clinical Laboratory, Approved Guideline, CLSI publication C28-A, Villanova, PA(1995).
5.Tietz, N. W., ed., Fundamentals of Clinical Chemistry, 3rd Edition, W.B. Saunders, Philadelphia, PA (1987).
6.CLSI. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP5-A2 [ISBN 1-56238-542-9. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087 USA, 2008.

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© 2015 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tüm hakları saklıdır
Üretici:Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
AT Temsilcisi:Shanghai International Holding Corp.GmbH(Avrupa)